

### **Endoscopic Applicator Indication**

The Endoscopic Applicator is intended for use in delivering hemostatic agent to bleeding surgical sites through a 5 mm or larger trocar.

### **Endoscopic Applicator Important Risk Information**

Do not use if damaged.

Do not use air as a propellant to extrude hemostatic agents instead of the stylet.

Do not use Endoscopic Applicator to manipulate or retract organs or tissue.

The Endoscopic Applicator should only be used by persons adequately trained in Minimal Invasive procedures.

Do not reuse or resterilize. Reuse of single- use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient.

Rx only: For safe and proper use of this device, please refer to full device Instructions for Use.