
FLOSEAL SPECIAL APPLICATOR TIPS INDICATION:

The FLOSEAL Malleable Tip and Trimmable Tip are intended for delivery of FLOSEAL Hemostatic Matrix.

Important Risk Information for FLOSEAL Special Applicator Tips

Do not use the FLOSEAL Malleable Tip and Trimmable Tip if the packaging is visibly damaged or if the tips are cracked.

Do not use air to remove residual FLOSEAL from the tip.

The FLOSEAL Malleable Tip and Trimmable Tip should only be used with FLOSEAL Hemostatic Matrix.

Read instructions for FLOSEAL Hemostatic Matrix prior to use.

The Malleable Tip should not be cut.

The Trimmable Tip should be trimmed away from the surgical cavity. Excess tip material should be collected in the inner pouch to discard.

Do not resterilize. Single-use only. Immediately discard after use.

RX only: For safe and proper use of these devices, please refer to full device Instructions for Use.